



Our new packaging concept

Many advantages
at the same price!

General Implants



- ✓ **Double sterile barrier**
Sterility for 5 years
- ✓ **Improved handling**
For caregivers and assistants
- ✓ **Direct product removal**
with screwdriver
- ✓ **Colour-coded labelling**
A handle for the correct diameter
- ✓ **Sealed outer packaging**
with practical tear tab

General Implants



General Implants

GmbH Deutschland

Max-Planck-Straße 39

78052 Villingen-Schwenningen

Germany


Phone +49 (7721) 99822-42

Fax +49 (7721) 99822-55


info@general-implants.com

www.general-implants.com

The procedure – step by step



① Open secondary packaging via tear tab and remove sterile blister



② Remove Tyvek foil from the outer blister and take out inner blister



③ Remove Tyvek foil from the inner blister



④ Insert screwdriver directly into the implant and remove



⑤ Take off protective cover for cover screw



⑥ Remove cover screw via screwdriver